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

ARCHITECT i2000_{SR}

- Fast Turnaround Time that delivers:
 - High Productivity
 - Up to 200 tests/hour
 - 135 sample load up (35stat 100routine)
 - 25 reagent positions refrigerated
- Confidence of Results:
 - **Chemiflex technology**
 - Pressure Diff - Clot detection
- Retest Sample Handler (RSH):
 - Maximizes operator efficiency
 - Multi-Dimensional Sampling
 - Automated Rerun, Dilution and Pretreatment
 - Easy sample access

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
Sample Carrier

- **Sample Carrier**
- **5 sample positions**
- **Allows different tubes of different sizes**
- **Primary or Aliquot tubes**
- **Sample Cups**

Sample Cup

Architect Reagents

Reagents

Calibrators

Controls



Processing Module: Reagent Center

- Reagent Kits (100 & 500 test kits size)
- Color coded , 2-D barcoded reagents
 - assay name
 - lot number
 - test size
 - expiration date
 - master cal curve info
- Septum is placed on reagent bottles prior to loading on the module
- Key Benefit: Easy inventory management
Easy to use.....simplify operations


Chemiluminescent Microparticle Immunoassay (CMIA)

A chemiluminescent label conjugated to the antibody or antigen, and it produces light when combined with its substrate. This method is very similar to MEIA, though the chemiluminescent reaction offers high sensitivity and ease of measurement. A noncompetitive sandwich format yields results that are directly proportional to the amount of analyte present.

ARCHITECT[®] uses a pre-trigger/trigger reaction to both enhance and generate light emission


Pre-Trigger is added


The Pre-Trigger is designed to release the Patented Acridinium away from the Solid Phase and into solution.


A magnet pulls the Paramagnetic Microparticles to the side of the Reaction Vessel wall — preventing the Microparticles from interfering in the light release/capture process.

Trigger is added to the solution — providing an alkaline medium for light generation.

The Acridinium label is released into solution.

The light generated is measured as **Relative Light Unit (RLU)**.

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
ARCHITECT HBsAg Assay Format (Qualitative)-1-step

ARCHITECT HBsAg Assay Format (Qualitative)

RESULTS

Calculations

- The ARCHITECT *i* System calculates the cutoff RLU from the mean RLU of three replicates of Calibrator 2. The cutoff RLU is stored for each reagent lot calibration.

$$\text{Cutoff RLU} = 0.013 \times (\text{Calibrator 2 mean RLU})$$

- The ARCHITECT *i* System calculates the S/CO result for each sample as follows:

$$\text{S/CO} = \text{Sample RLU/Cutoff RLU}$$

Interpretation of Results

ARCHITECT HBsAg Initial Results	
Instrument Interpretation	Retest Procedure
Nonreactive (S/CO value < 1.00)	No retest required.
Reactive (S/CO value \geq 1.00)	Retest in duplicate.

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
S/CO ≥ 1.0 = Reactive

ARCHITECT HCV-Ab Assay Format


HCr43: The HCr43 protein is expressed in *Escherichia coli* (*E. coli*) and is composed of two noncontiguous coding regions of the HCV genome sequence. The first region represents amino acids 1192 to 1457 (33c) of the HCV sequence. The second of the two regions represents amino acids 1 to 150 (core) of the HCV sequence. Because of the similarity of the genomic organization of the flaviviruses, it is suggested that the first sequence is from the NS3 coding region and the second sequence is from the core coding region of HCV.

c100-3: The c100-3 antigen is a recombinant HCV protein expressed in *Saccharomyces cerevisiae* (yeast). The genomic organization of flaviviruses suggests that the cloned sequence is contained within the putative nonstructural (NS3 and NS4) regions of HCV. The c100-3 protein is a chimeric fusion protein with 154 amino acids of human superoxide dismutase (hSOD), five linker amino acids, amino acids number 1569 to 1931 of the HCV polyprotein, and the additional five amino acid linker at the carboxyl terminus.


ARCHITECT HCV-Ab Assay Format

Assay Schematic


ARCHITECT HCV-Ab Assay Format

SAMPLE CAPTURE

SAMPLE DETECTION

ARCHITECT HCV-Ab Assay Format

RESULTS

Calculation

- The ARCHITECT *i* System calculates the cutoff RLU from the mean chemiluminescent signal of three Anti-HCV Calibrator 1 replicates and stores the result. The cutoff RLU is determined by multiplying the Anti-HCV Calibrator 1 mean RLU by 0.074.

$$\text{Cutoff RLU} = \text{Calibrator 1 Mean RLU} \times 0.074$$

- The ARCHITECT *i* System calculates a result based on the ratio of the sample RLU to the cutoff RLU (S/CO) for each specimen and control.

$$\text{S/CO} = \text{Sample RLU/Cutoff RLU ratio}$$

Interpretation of Results

Initial ARCHITECT Anti-HCV Results			
Initial Result (S/CO)	Instrument Flag	Interpretation	Retest Procedure
≥ 1.00	REACTIVE	Reactive	No retest required.
0.80 to 0.99	GRAYZONE	Grayzone	Retest in duplicate.
0.00 to 0.79	NONREACTIVE	Nonreactive	No retest required.

ARCHITECT HIV Ag/Ab Combo Assay Format (1)

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S/CO < 1.0 = Nonreactive


S/CO ≥ 1.0 = Reactive

ARCHITECT HIV Ag/Ab Combo Assay Format (2)

Anti-HIV Antibody Detection

Microparticle

rHIV-1 Group M Envelope Protein
rHIV-2 Envelope Protein
rHIV-1 Group O Envelope Protein

Conjugate

rHIV-1 Group M Envelope Protein Conjugate

HIV-2 Peptide Conjugate


HIV-1 Group M Peptide Conjugate

rHIV-1 Group O Envelope Protein Conjugate

ARCHITECT HIV Ag/Ab Combo Assay Format (3)

HIV p24 Antigen Detection

Microparticle
Monoclonal Ab
(Conformational)

Conjugate
Monoclonal Ab
(Conformational)

100 μ L Sample
50 μ L Specimen Diluent
50 μ L Microparticle

18 min
Incubation

Wash

Add 50 μ L Conjugate


4 min
Incubation

Wash

Pretrigger
Trigger
Report Result

ARCHITECT HIV Ag/Ab Combo Assay Format (4)

Note: The two types of HIV coated microparticles are blended during manufacturing and filled into the microparticle reagent bottle. The five conjugate types are also blended during manufacturing and filled into the conjugate reagent bottle.

ARCHITECT Syphilis Ab Assay Format

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S/CO < 1.0 = Nonreactive

S/CO ≥ 1.0 = Reactive

ARCHITECT Syphilis Ab Assay Format

Calculation

The ARCHITECT Syphilis TP assay calculates a result based on a cutoff determined by the following calculation.


- $\text{Cutoff (CO)} = \text{Calibrator 1 Mean RLU} \times 0.20$
- $\text{S/CO} = \text{Sample RLU} / \text{Cutoff RLU}$
- The cutoff RLU is stored for each reagent lot calibration.

Interpretation of Results

- Specimens with S/CO values < 1.0 are considered nonreactive by the ARCHITECT Syphilis TP assay.
- Specimens with S/CO values > 1.0 are considered reactive by the ARCHITECT Syphilis TP assay.


ARCHITECT Syphilis Ab Assay Format

Assay Schematic

ARCHITECT Syphilis Ab Assay Format

ARCHITECT Syphilis TP Assay Format

The End